



Leica EM IGL

Automated Immunolabelling System
according to Posthuma

Immunolabelling, tedious in the past ...

With the introduction of particulate immuno-markers in the 1970's and its revolution in 1978 when Roth used colloidal gold particles as a marker, a new era began in the localization of the constituents of the cell. Immunogold labelling became the method of choice for detection of antigens in biological samples.

Until today, immunogold labelling of ultrathin sections on grids has been carried out completely by hand.

- Drops of different reagents are arranged in a highly specific order on a very clean surface.
- Each grid is carefully transferred from one drop to the other – by hand.
- The sequence of grid transfer must be controlled very carefully to avoid mistakes.
- The duration of each immunoreagent and wash step must be carefully noted to achieve reproducible results.

But why not reduce your effort?

- Why spend many hours in front of a Parafilm sheet?
- Why risk losing grids/sections while transferring them from one drop to the other?
- Why risk ruining a complete run by transferring a grid to the wrong droplet?
- Why risk cross contamination by contaminated forceps or loop?
- Why have results with low reproducibility?
- Why not turn this into an automatic process?

... routine today!



LEICA EM IGL Designed by Werner Hölbl

The LEICA EM IGL meets the demands of today's immunolabelling laboratory. By preparing all your reagents together, you can set up a run within a few minutes and free your time for more productive activities. Reagent droplets are placed in humidified carries and moved automatically to grids according to the programme in use.

The Leica EM IGL gives you

- Labelling of 24 grids simultaneously.
- The ability to use up to 24 different primary antibodies in one run.
- 80% time saving compared to the manual method.
- Cost saving - you are free for more productive activities.
- Specimen safety – the correct sequence of reagents will always be applied.
- Minimising cross contamination from forceps and loops during labelling.
- Reproducible results from each immunolabelling run by following the exact incubation times according to your protocol.
- The possibility of standard contrasting with lead citrate and uranyl acetate.

From easy programming ...


Latest technology and PC programming software is used to interface user and instrument.
A simple to use control panel allows data storage of

- Up to 99 programmes.
- All of which can be named.
- Reagents can be selected from a reagent list which can be customised by the user.

Programming is carried out either via the membrane-covered keyboard or PC.
Audible warnings and clear hints are given to allow simple programming.
Intelligent software allows labelling to continue after a power failure.



1



Display

(The example shows the screen during programming).

- Programme number and name
- Slide number and name
- Duration time, droplet size and continue (CONT) or PAUSE before the next slide
- Total time of programme

2



Selection and adjusting of programmes

3



Keys to start and pause a programme

4



Selecting a programme
Interrupting a programme

5



Manual movement of a slide
Keyboard lock-function

... via easy set up ...

Load up to 24 grids onto the magnetic holder
Load droplets onto slides
Load into EM IGL
Press START

Logical droplet management pads provide clear control over the droplet and grid arrangement, thus preventing any user error. A special Teflon-coated grid holder accepts up to 24 nickel grids and ensures the back of the grid remains absolutely dry.



Date	02.09.2003	
Operator		
Experiment no		
Program no		
Program name		

Gridholder	Pos.	Specimen
	1 a	MDCK Cells
	1 b	MDCK Cells
	1 c	MDCK Cells
2	2 a	MDCK Cells
	2 b	MDCK Cells
	2 c	MDCK Cells
3	3 a	MDCK Cells
	3 b	MDCK Cells
	3 c	MDCK Cells
4	4 a	Rat pancreas
	4 b	Rat pancreas
	4 c	Rat pancreas
5	5 a	Rat pancreas
	5 b	Rat pancreas
	5 c	Rat pancreas
6	6 a	Rat pancreas
	6 b	Rat pancreas
	6 c	Rat pancreas
7	7 a	Rat liver
	7 b	Rat liver
	7 c	Rat liver
8	8 a	Rat liver
	8 b	Rat liver
	8 c	Rat liver

Loading the grids onto the holder in the correct position according to the droplet management pad.

Special '24 well' slides hold the reagent drops in exactly the correct position in relation to the grids. Immunoreagents are 5 µl in volume and rinse steps are 30 µl. For quick addition of wash and buffer droplets to the slides, a multi 8-channel pipette is provided. A single channel pipette is used for adding the valuable antibody drops.



Using a multi 8-channel pipette for fast and accurate addition of wash and buffer droplets.

Leica EM IGL
Droplet Management Pad

Slide	Pos.	Immuno Reagent 1	Immuno Reagent 2
1	a	MAB CD 63	
	b	MAB CD 63	Ratcbl anti mouse
	c	MAB CD 63	Ratcbl anti mouse
2	a	Rx M MHC II	
	b	Rx M MHC II	Ratcbl anti mouse
	c	Rx M MHC II	
3	a	RxM MHC I	
	b	RxM MHC I	
	c	RxM MHC I	
4	a	RxRat Amylase	
	b	RxRat Amylase	
	c		
5	a	Rab T7	
	b	Rab T8	Ratcbl anti mouse
	c	Rab T9	Ratcbl anti mouse
6	a	Rab 20	Ratcbl anti mouse
	b	Rab 21	Ratcbl anti mouse
	c	Rab 22	Ratcbl anti mouse
7	a	MAB A50P 1/60	Ratcbl anti mouse
	b	MAB A50P 1/60	Ratcbl anti mouse
	c		Ratcbl anti mouse
8	a	MFR 46 1/50000	
	b	MFR 46 1/50000	Ratcbl anti mouse
	c	MFR 46 1/50000	



Using a single channel pipette to add the antibody drops.



... and easy loading ...



Loaded grid holder with grids facing down is placed into the grid station of the LEICA EM IGL.

Beginning with the first reagent according to your programme the stack of the LEICA EM IGL is loaded with the slide carriers.



Once the slides and grids are loaded to the LEICA EM IGL and START button is pressed **YOU can carry out more important things than moving grids from one drop to the next.**

... to perfect results on cryo sections ...

Rat pancreas.


Bar = 200 nm

S = Secretory granule

G = Golgi complex

R = Rough endoplasmatic reticulum

Courtesy of: George Posthuma
University Medical Centre Utrecht, Netherlands



HELA cell.


Incubated with cholesterol
marker coupled to biotin.

Bar = 200 nm

P = Plasma membrane

E = Endosomes/lysosome

Courtesy of: George Posthuma
University Medical Centre Utrecht, Netherlands



Rat liver,


peroxisomes labelled with rabbit anti-catalase antibodies.

Detected with goat anti-rabbit ultra small gold particles.


Silver enhanced with R-Gent SE-EM.

Bar = 1 μ m

Courtesy of: Bruno M. Humbel
Molecular Cell Biology
Faculty of Biology,
Utrecht University,
Netherlands



... and perfect results on resin sections ...




Mouse pancreas

Cryo fixed followed by freeze substitution and low temperature embedded in HM20. Immunogold/silver labelling of alpha-amylase using Aurion Ultrasmall gold conjugate and R-gent SE-EM silver enhancement reagent.

Bar = 2 μ m

Courtesy of: Hong Yi, Emory School of Medicine
Microscopy Core, Emory Univ., Atlanta, USA.

... to viruses on coated grids.



Sucrose gradient purified Modified Virus Ankara (MVA).

The B5R envelope protein is labelled with a rabbit anti-B5R antibody and detected with a goat anti-rabbit antibody coupled to 6 nm gold particles.

Bar = 200 nm

Courtesy of: Daniele Spehner
INSERM, Strasbourg, France.

Technical Specifications:

Dimensions:

Width: 50 cm
Depth: 26 cm
Height: 54 cm
Weight: 16 kg

Electrical power consumption: 50 W
Wide range power supply for 100 – 240 VAC, 50 – 60 Hz.
Power is via standard socket. No special room requirements are needed.

Set up:

With the order number 70 89 01 a complete working outfit with all necessary accessories and consumables will be delivered (except reagents for immunogold labelling and cleaning).
For the cleaning of glass slides we recommend 1 % DECON 90.

Order number 70 89 01 consists of:

LEICA EM IGL Basic instrument.
Slide carrier blue (30 pcs.)
Slide carrier brown (30 pcs.)
Filter paper (300 pcs.)
Glass slides (100 pcs.)
Grid holder (1 pc.)
Nickel grids (100 pcs.)
Nickel grids carbon/formvar coated, 100 mesh hexagonal (50 pcs.)
Pipette multi 8-channel, 10–100 µl, 1 pc. (96 tips included).
Pipette single channel, 0.5-10 µl, 1 pc. (96 tips included)
Fine forceps with straight tips, non magnetic (1 pc.)
Plastic forceps (1 pc.)
Reagent tub and slide rack for cleaning glass slides
Droplet management pad
Programme protocol pad
Leica EM IGL – PC programming software package

Leica Microsystems – the brand for outstanding products

Leica Microsystems' mission is to be the world's first-choice provider of innovative solutions to our customers' needs for vision, measurement, lithography and analysis of microstructures.

Leica, the leading brand for microscopes and scientific instruments, developed from five brand names, all with a long tradition: Wild, Leitz, Reichert, Jung and Cambridge Instruments. Yet Leica symbolizes innovation as well as tradition.

Leica Microsystems – an international company with a strong network of customer services

Australia:	Gladesville	Tel. +61 2 9879 9700	Fax +61 2 9817 8358
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison	Tel. +33 1 473 285 85	Fax +33 1 473 285 86
Germany:	Bensheim	Tel. +49 6251 136 0	Fax +49 6251 136 155
Italy:	Milan	Tel. +39 0257 486.1	Fax +39 0257 40 3273
Japan:	Tokyo	Tel. +81 3 5435 9600	Fax +81 3 5435 9615
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Sollentuna	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Glattbrugg	Tel. +41 1 809 34 34	Fax +41 1 809 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives of Leica Microsystems
in more than 100 countries.

The companies of the Leica Microsystems Group operate internationally in five business segments, where we rank with the market leaders.

• Microscopy Systems

Our expertise in microscopy is the basis for all our solutions for visualization, measurement and analysis of microstructures in life sciences and industry. With confocal laser technology and image analysis systems, we provide three-dimensional viewing facilities and offer new solutions for cytogenetics, pathology and materials sciences.

• Specimen Preparation


We provide comprehensive systems and services for clinical histo- and cytopathology applications, biomedical research and industrial quality assurance. Our product range includes instruments, systems and consumables for tissue infiltration and embedding, microtomes and cryostats as well as automated stainers and coverslippers.

• Medical Equipment

Innovative technologies in our surgical microscopes offer new therapeutic approaches in microsurgery. With automated instruments for ophthalmology, we enable new diagnostic methods to be applied.

• Semiconductor Equipment

Our automated, leading-edge measurement and inspection systems and our E-beam lithography systems make us the first choice supplier for semiconductor manufacturers all over the world.



www.em-preparation.com

Leica Mikrosysteme GmbH Telephone +43 1 48899
Hernalser Hauptstrasse 219 Fax +43 1 48899-350
A-1170 Vienna, Austria

Leica
MICROSYSTEMS