

此表中质粒均为平台所有，可提供常用工具AAV包装

No.	Catalog	Item Name
1	20296	pAAV-EF1a-double floxed-EYFP-WPRE-HGHpA
2	20297	pAAV-EF1a-double floxed-hChr2 (H134R) -mCherry-WPRE-HGHpA
3	26966	pAAV-Ef1a-DIO eNpHR 3.0-EYFP
4	26969	pAAV-CaMKIIa-hChr2 (H134R) -EYFP
5	26971	pAAV-CaMKIIa-eNpHR 3.0-EYFP
6	26972	pAAV-hSyn-eNpHR 3.0-EYFP
7	26976	pAAV-hSyn-hChr2 (H134R) -mCherry
8	28305	pAAV-FLEX-ArchT-tdTomato
9	35498	pAAV-Ef1a-DIO C1V1 (t/t)-TS-mCherry
10	41015	hSyn1-FLEX-hChr2-tdtomato
11	44361	pAAV-hSyn-DIO-hM3D (Gq) -mCherry
12	44362	pAAV-hSyn-DIO-hM4D (Gi) -mCherry
13	50459	pAAV-hSyn-DIO-mCherry
14	50460	pAAV-EF1a-DIO-hM3D (Gq) -mCherry
15	50461	pAAV-EF1a-DIO-hM4D (Gi) -mCherry
16	50462	pAAV-EF1a-DIO-mCherry
17	50465	pAAV-hSyn-EGFP
18	50469	pAAV-CaMKIIa-EGFP
19	50473	pAAV-GFAP-EGFP
20	50474	pAAV-hSyn-HA-hM3D (Gq) -mCherry
21	50475	pAAV-hSyn-HA-hM4D (Gi) -mCherry
22	50476	pAAV-CaMKIIa-hM3D (Gq) -mCherry
23	50477	pAAV-CaMKIIa-hM4D (Gi) -mCherry
24	50942	pAAV-hSyn1-mRuby2-GSG-P2A-GCaMP6s-WPRE-pA
25	50943	pAAV-hSyn1-mRuby2-GSG-P2A-GCaMP6f-WPRE-pA
26	50955	AAV-flex-ReaChR-citrine
27	51082	AAV-EF1a-DIO-GCaMP6s-P2A-nls-dTomato
28	51083	AAV-EF1a-DIO-GCaMP6f-P2A-nls-dTomato
29	51084	AAV-hSyn1-GCaMP6s-P2A-nls-dTomato
30	51085	AAV-hSyn1-GCaMP6f-P2A-nls-dTomato

31	51086	AAV-CaMKIIa-GCaMP6s-P2A-nls-dTomato
32	51087	AAV-CaMKIIa-GCaMP6f-P2A-nls-dTomato
33	51473	pAAV-hSyn1-mRuby2-GSG-P2A-GCaMP6m-WPRE-pA
34	55632	pAAV-Efla-mCherry-IRES-Cre
35	55636	pAAV-EF1a-Cre
36	55637	pAAV-EF1a-Flpo
37	59462	pAAV-CAG-tdTomato (codon diversified)
38	63141	pAAV-Syn-ArchT-GFP
39	83279	pscAAV-CAG-GFP
40	EP124	CamKII-mcherry-cre
41	EP126	synapsin-mcherry-cre
42	X888	hSyn-CRE-EGFP
43	GP73	pAAV-hsyn-DIO-GCamp6s
44	GP74	pAAV-hSyn-GCaMP6s
45	GP76	pAAV-CAG-EGFP
46	GP77	pscAAV-CAG-Cre
47	GP78	pscAAV-CAG-mCherry
48	GP79	pAAV-EF1a-DIO-EGFP
49	GP80	pAAV-hsyn-DIO-GCamp6f
50	GP81	pAAV-hsyn-GCamp6f
51	GP82	pAAV-hsyn-Cre
52	GP84	pAAV-hsyn-mcherry
53	GP85	XP331_pX601-AAV-EFs-Cas9-SV40-ployA(+taa)
54	GP86	pAAV hSyn-CRE-P2A-mcherry
55	GP87	GRP024-GFAP-tdTomato-WPRE
56	GP91	pCX1519_Lenti_EF1 α -EGFP-P2A-Puro-invU6-BpiI-BGHpA